

UMBRA

Premium Tattoo & Body Piercing Nitrile Examination Gloves



- Latex free
- Powder free
- Disposable
- Non-Sterile
- Ambidextrous
- Beaded Cuff



BASIC MATERIALS

Synthetic nitrile rubber free of any latex protein

DONNING POWDER

Free of powder (free of talcum, magnesium silicate) according to EN 455-3

FEATURES

- Latex free
- Powder free
- Disposable
- Non-Sterile
- Ambidextrous
- Beaded Cuff

INTENDED USE

- Tattooing
- Piercing
- Patient Examination
- Food Service
- Hygiene
- Beauty & Salons

PHYSICAL DIMENSION AND WEIGHT

SIZE	LENGTH (MM)	PALM WIDTH (MM)	THICKNESS SINGLE WALL (MM)		
			FINGER	PALM	CUFF
X-Small	240 +/-10	80 +/-5	0.12 +/-0.03	0.08 +/-0.03	0.07 +/-0.03
Small	240 +/-10	85 +/-5	0.12 +/-0.03	0.08 +/-0.03	0.07 +/-0.03
Medium	240 +/-10	95 +/-5	0.12 +/-0.03	0.08 +/-0.03	0.07 +/-0.03
Large	240 +/-10	105 +/-5	0.12 +/-0.03	0.08 +/-0.03	0.07 +/-0.03
X-Large	240 +/-10	115 +/-5	0.12 +/-0.03	0.08 +/-0.03	0.07 +/-0.03

PHYSICAL PROPERTIES

CHARACTERISTICS	REFERENCE STANDARDS	BEFORE AGING	AFTER AGING
Force at Break (N)	EN 455-2	Min 4 N	Min 3.8 N
Tensile Strength (MPa)	ASTM: D6319	Min 14 MPa	Min 14 MPa
Elongation at Break (%)	ASTM: D6319	Min 500%	Min 400%

QUALITY ASSURANCE

CHARACTERISTICS	INSPECTION LEVEL & SAMPLING PLAN	AQL & ACCEPTANCE LIMIT
Major Defects	S-4	AQL 1.5
Minor Defects	S-4	AQL 1.5
Physical Dimension	S-4	AQL 1.5

PACKING AND STORAGE INSTRUCTIONS

Packing	Typically 100 gloves per inner box; 10 boxes per outer carton
Humidity Range	Keep in a cool dry place
Temperature Range	Store at temperature no less than 5°C and no more than 40°C
Sunlight Effect	Keep away from direct sunlight or fluorescent light
Shelf Life	ASTM D7160 Compliant. Shelf life is 5 years from the manufacture date under proper storage conditions.

RECOVERYAFTERCARE.COM

Distributed by Recovery Inc.
© Recovery, Inc. 2019

Questions or Comments?
support@recoveryaftercare.com